



TripleNine Stainless is the distributor of ITT Pure-Flo products in NSW, Australia. ITT Pure-Flo offers superior precision-engineered hygienic diaphragm valves and process components for the pharmaceutical, bioprocessing and fine chemical industries.

Through both standard and custom designed valve assemblies, ITT Pure-Flo is committed to providing the best quality and value in engineered solutions for your unique flow-control needs

As a recognized leader in innovation for the Pharmaceutical and Bioprocessing Industries, ITT prides it self in developing solutions for process problems. Accomplished through direct client input, this has, and will continue to be, our focus.

Please contact us for advice on ITT valves that are best suited for your application.



## Valve Types



[Two Way Diaphragm Valve](#)



[Bio-Tek Valve](#)



[Zero Static Block Body](#)



[Sterile Access and GMP Valve](#)



[Integral Sterile Access Valves](#)



[Zero Static Back to Back Sample Valve](#)



[Multiport Divert Valves](#)



[Tank Bottom Valve](#)



[Sample and Bleed Valve](#)

•



[Zero Static with Upstream Sample & Downstream Purge](#)

•



[Dual Flow and Bypass](#)

•



[Zero Static Down Stream Purge](#)


•



[Chromatography Valve](#)



[Custom Block Designs](#)



[Modular Valve Assemblies](#)



[Sterile Barrier, Block and Bleed Valve](#)



[Sterile Filter Shunt Valve](#)



[Specialty Zero Static Use Points](#)

## **Service and Parts**

Ensuring that process valves are kept in good operating condition is critical in the highly regulated Pharmaceutical industry. Routine preventative maintenance will pay many returns in uptime and trouble free operation.

TripleNine Stainless is committed to supporting customers with valve service and the supply of quality authentic replacement parts such as Diaphragms, critical spare topworks, gall resistant fasteners, and autoclave rebuild kits are fundamental components of the service parts program.

All ITT service parts meet our strict quality standards. Authentic ITT service parts are identical to those used in our valve production, giving customers peace of mind that the part they receive will be equal to the part they are replacing. Please contact us for more details on service and parts.

